

cobas p 630 Instrument

The pre-analytics solution that makes life easier



The **cobas p 630** Instrument offers in combination with the COBAS® AmpliPrep/ COBAS® TaqMan® System a fully automated pre-analytical solution for primary tube handling. The system automatically pipettes primary and secondary tubes and controls into sample input tubes for the COBAS AmpliPrep Instrument.

The **cobas p 630** Instrument can be combined with up to three COBAS AmpliPrep Instruments and AmpliLink® software to ensure full traceability of workflow.



cobas p 630 Instrument

Your benefit

Efficiency

- Automated handling of primary and secondary tubes

Flexibility

- Compatible with a variety of sample tubes
- Modular design

Full traceability

- Barcode tracking from patient tube to result

Process surveillance

- Monitors liquid handling

Product characteristics

- Uncapping and recapping of the sample tube
- Pipetting Roche controls from control tubes to sample tubes
- Pipetting samples from primary and secondary tubes to sample tubes
- Multiple tests can be ordered on a single primary tube
- Only one LIS interface required

Unit dimensions

- 112 cm wide, 101 cm deep, 90 cm high

Sample processing throughput

- 320 samples on board
- 154 tubes per hour for 650 µL samples
- 149 tubes per hour for 1.0 mL samples
- 157 tubes per hour for 500 uL samples

